

Versa**TREK**[®]

Media Specifications

Blood Culture

TREK Media*	Catalog Number	Specimen Type	Specimen Volume
VersaTREK REDOX 1® 80 ml with stir bar	7102-44	Blood Sterile Body Fluid	Range 0.1 to 10 ml
VersaTREK REDOX 2® 80 ml	7103-44	Blood Sterile Body Fluid	Range 0.1 to 10 ml
VersaTREK REDOX 1 EZ Draw® 40 ml with stir bar	7106-44	Blood Sterile Body Fluid	Range 0.1 to 5 ml
VersaTREK REDOX 2 EZ Draw 40 ml	7107-44	Blood Sterile Body Fluid	Range 0.1 to 5 ml
VersaTREK REDOX 2M EZ Draw 40 ml	7110-44	Blood Sterile Body Fluid	Range 0.1 to 5 ml

*50 bottles/case



Mycobacteria

TREK Media	Catalog Number	Specimen/Product Description
VersaTREK Myco**	7111-42	Blood Sterile Body Fluid Processed Specimens
VersaTREK Myco GS	7112-42	Myco Growth Supplement
VersaTREK Myco PVNA	7113-42	Myco Antibiotic Supplement
VersaTREK Myco AS	7114-42	Myco Antibiotic Supplement
VersaTREK Myco Susceptibility Kit (80 tests/kit)	7115-60	Rifampin, Isoniazid, Ethambutol
Myco PZA Kit (25 tests/kit)	7116-70	Pyrazinamide



**Specimen volume- up to 1 ml

VersaTREK[®]

240 Model



Dimensions

52 in. width, 30 3/4 in. depth
40 1/4 in. height
132 cm width, 78 cm depth
103 cm height

Unit Weight (Empty)

656 lbs
304 kg

Unit Weight (Fully Loaded)

766 lbs
348 kg

Capacity

96 to 240 bottles

528 Model



Dimensions

52 in. width, 30 3/4 in. depth
76 1/4 in. height
132 cm width, 78 cm depth
194 cm height

Unit Weight (Empty)

1,212 lbs
561 kg

Unit Weight (Fully Loaded)

1,455 lbs
660 kg

Capacity

Up to 528 bottles

VersaTREK Windows[®] Software

VersaTREK Computer Package* 6110-31

* The VersaTREK Computer Package includes 15" Touch Screen Monitor, Microsoft Access XP, VersaTREK Windows Software CD, HP Laserjet Printer, UPS 1200 VA Universal 6-Outlet, and VersaTREK Operator's CD Manual.

VersaTREK Accessories

VersaTREK Myco Adaptors, 25/box.....	6040-30	Thermometer, 18°C - 50°C, 135mm, each*	6440-30
Red Disable Location Cap, each*.....	6044-30	Absorbent Sheet, each*.....	6444-30
VersaTREK Stylus, 4/pack*.....	6420-30	Label, Drawer, Numbered 1 - 22, each*.....	6452-30
VersaTREK Stylus Holder, each.....	6422-30	VersaTREK EZ-ID Labels, 100/pack.....	6450-35
Bottle Carry Rack, 24-positions*.....	6424-30	VersaTREK Connectors, 50/box.....	7150-44

* These items are included in the VersaTREK Starter Kit (6499-31). Quantity of items in the Starter Kit are: 5 red disable location caps, 1 VersaTREK stylus 4-pack, 1 bottle carry rack, thermometer, 3 absorbent sheets, 1 set of drawer labels. One cable, CAT5E 25' Snagless Patch (6172-30), is also included with the Starter Kit. One complete Starter Kit is shipped with each VersaTREK instrument.

Packaging and Storage Information

- All VersaTREK REDOX 80 ml and VersaTREK REDOX EZ Draw 40 ml bottles are packaged 50/box and should be stored at 15-30°C away from light.
- VersaTREK Myco contains 50/box and should be stored at 15-30°C, away from light.
- VersaTREK Myco GS, Myco AS and Myco PVNA are each packaged 5/box and should be stored 2-8°C. Myco GS must be stored in the dark.
- For technical information, refer to the technical inserts: VersaTREK REDOX, VersaTREK Myco, VersaTREK Myco Susceptibility and VersaTREK PZA.