

# Versa**TREK**<sup>®</sup>

## Media Specifications

### Blood Culture

TREK Media*	Catalog Number	Specimen Type	Specimen Volume
VersaTREK REDOX 1® 80 ml with stir bar	7102-44	Blood Sterile Body Fluid	Range 0.1 to 10 ml
VersaTREK REDOX 2® 80 ml	7103-44	Blood Sterile Body Fluid	Range 0.1 to 10 ml
VersaTREK REDOX 1 EZ Draw® 40 ml with stir bar	7106-44	Blood Sterile Body Fluid	Range 0.1 to 5 ml
VersaTREK REDOX 2 EZ Draw 40 ml	7107-44	Blood Sterile Body Fluid	Range 0.1 to 5 ml
VersaTREK REDOX 2M EZ Draw 40 ml	7110-44	Blood Sterile Body Fluid	Range 0.1 to 5 ml

\*50 bottles/case



### Mycobacteria

TREK Media	Catalog Number	Specimen/Product Description
VersaTREK Myco**	7111-42	Blood Sterile Body Fluid Processed Specimens
VersaTREK Myco GS	7112-42	Myco Growth Supplement
VersaTREK Myco PVNA	7113-42	Myco Antibiotic Supplement
VersaTREK Myco AS	7114-42	Myco Antibiotic Supplement
VersaTREK Myco Susceptibility Kit (80 tests/kit)	7115-60	Rifampin, Isoniazid, Ethambutol
Myco PZA Kit (25 tests/kit)	7116-70	Pyrazinamide



\*\*Specimen volume- up to 1 ml

# VersaTREK<sup>®</sup>

## 240 Model



### Dimensions

52 in. width, 30 3/4 in. depth  
40 1/4 in. height  
132 cm width, 78 cm depth  
103 cm height

### Unit Weight (Empty)

656 lbs  
304 kg

### Unit Weight (Fully Loaded)

766 lbs  
348 kg

### Capacity

96 to 240 bottles

## 528 Model



### Dimensions

52 in. width, 30 3/4 in. depth  
76 1/4 in. height  
132 cm width, 78 cm depth  
194 cm height

### Unit Weight (Empty)

1,212 lbs  
561 kg

### Unit Weight (Fully Loaded)

1,455 lbs  
660 kg

### Capacity

Up to 528 bottles

## VersaTREK Windows<sup>®</sup> Software

VersaTREK Computer Package\* ..... 6110-31

\* The VersaTREK Computer Package includes 15" Touch Screen Monitor, Microsoft Access XP, VersaTREK Windows Software CD, HP Laserjet Printer, UPS 1200 VA Universal 6-Outlet, and VersaTREK Operator's CD Manual.

## VersaTREK Accessories

VersaTREK Myco Adaptors, 25/box.....	6040-30	Thermometer, 18°C - 50°C, 135mm, each* .....	6440-30
Red Disable Location Cap, each*.....	6044-30	Absorbent Sheet, each*.....	6444-30
VersaTREK Stylus, 4/pack*.....	6420-30	Label, Drawer, Numbered 1 - 22, each*.....	6452-30
VersaTREK Stylus Holder, each.....	6422-30	VersaTREK EZ-ID Labels, 100/pack.....	6450-35
Bottle Carry Rack, 24-positions*.....	6424-30	VersaTREK Connectors, 50/box.....	7150-44

\* These items are included in the VersaTREK Starter Kit (6499-31). Quantity of items in the Starter Kit are: 5 red disable location caps, 1 VersaTREK stylus 4-pack, 1 bottle carry rack, thermometer, 3 absorbent sheets, 1 set of drawer labels. One cable, CAT5E 25' Snagless Patch (6172-30), is also included with the Starter Kit. One complete Starter Kit is shipped with each VersaTREK instrument.

## Packaging and Storage Information

- All VersaTREK REDOX 80 ml and VersaTREK REDOX EZ Draw 40 ml bottles are packaged 50/box and should be stored at 15-30°C away from light.
- VersaTREK Myco contains 50/box and should be stored at 15-30°C, away from light.
- VersaTREK Myco GS, Myco AS and Myco PVNA are each packaged 5/box and should be stored 2-8°C. Myco GS must be stored in the dark.
- For technical information, refer to the technical inserts: VersaTREK REDOX, VersaTREK Myco, VersaTREK Myco Susceptibility and VersaTREK PZA.